

HOW TO USE THE AEROCHAMBER® WITH A MASK

1. Shake the metered dose inhaler (MDI) and remove its cap.
2. Insert the mouthpiece of the MDI into the end of the AeroChamber®.
3. Apply the mask to the child's face so there are no leaks between the child's face and the mask.
4. Have the child breathe out.
5. Push down once on the MDI canister.
6. Have the child inhale slowly and deeply and hold breath for 10 seconds or as long as comfortable. If this is not possible, have the child breathe normally from the device 5-6 times.
7. Replace the cap.

If a second dose is prescribed, repeat above steps, waiting one minute between each dose.

