

HOW TO USE THE AEROCHAMBER®

1. Shake the metered dose inhaler (MDI) and remove its cap.
2. Remove the AeroChamber® cap and insert the mouthpiece of the MDI into the end of the AeroChamber®.
3. Place the mouthpiece of the AeroChamber® in the child's mouth making sure there is a tight seal.
4. Have the child breathe out.
5. Push down once on the MDI canister.
6. Have the child inhale slowly and deeply and hold breath for 10 seconds or as long as comfortable. If this is not possible, have the child breathe normally from the device 5-6 times. It is recommended that the child holds breath, if possible.
7. Replace the caps.

If a second dose is prescribed, repeat above steps, waiting one minute between each dose.

